

### Introduction

This application note describes the general results when measuring animal feed using the Phoenix NIR Analyzers.

### Method

All samples were run on a Phoenix Series NIR Analyzer and calibrations were developed using Alligator Calibration software.

All samples were ground and mixed well before being analyzed.

The large rotating cup was used for analysis. Sample cup was ½ filled with lid.

### Results

Samples were collected from many locations around the world and include many different types of broiler formulas.

All results are on as-is basis.

### Definitions

**# of Samples:** Total number of samples for each constituent.

**Range:** Constituent range in the calibration

**R<sup>2</sup>:** Correlation coefficient is the agreement between the wet chemistry and NIR results. Correlation is dependent on lab accuracy and constituent range.

**SECV:** Cross Validation Error of the calibration. This value is approximately what can be expected when using the calibration for routine analysis.

Constituent	# of Samples	Range	R <sup>2</sup>	SECV
Protein	300	16.25 – 22.9	.919	.47
Moisture	301	9 – 16	.948	.32
Fat	300	2.2 – 8.5	.970	.27

