

### Introduction

This application note describes the general results when measuring cocoa butter using for the Phoenix NIR Analyzers.

### Method

All samples were run on a Phoenix Series NIR Analyzer and calibrations were developed using Alligator Calibration software.

All samples were analyzed between 40-50° C. Samples were mixed well before loading in the cup taking care to minimize air bubbles.

The following sample accessories were used:

Ring Cup

0.3 mm reflector

### Results

Samples were collected from many locations around the world and include many different types of cocoa butter.

All results are on as-is basis.

### Definitions

**# of Samples:** Total number of samples for each constituent.

**Range:** Constituent range in the calibration

**R<sup>2</sup>:** Correlation coefficient is the agreement between the wet chemistry and NIR results. Correlation is dependent on lab accuracy and constituent range.

**SECV:** Cross Validation Error of the calibration. This value is approximately what can be expected when using the calibration for routine analysis.

Constituent	# of Samples	Range	R <sup>2</sup>	SECV
FFA	472	0.1 – 2.2	0.97	0.077
Iodine Value (IV)	204	33.5 – 42.3	0.99	0.3

